### PMTACS – NY AMS METADATA REPORT

Frank Drewnick
Atmospheric Sciences Research Center, University at Albany
05/14/2002

### **Instrument Description**

The Aerosol Mass Spectrometer (AMS, Aerodyne Research Inc.) is a compact on-line measuring quadrupole mass spectrometer, based on thermal evaporation of aerosol particles.

The particles are focused by an aerodynamic lens while they are accelerated into a differentially pumped vacuum chamber. The inlet flow is 0.1 l/min. Particles in a diameter range from 40 nm to 600 nm are transmitted through the lens with 100 % efficiency. The particle beam passes a skimmer and a chopper before it crosses a flight chamber (0.39 m) and impacts onto a heated surface (600 – 900 °C). Volatile and semi volatile particle components are evaporated on the heater. The vapor is ionized by electron impact ionization and analyzed by a quadrupole mass spectrometer. Ions are detected with a multiplier. A more detailed description of the instrument is given in Jayne et al. (2000).

The instrument can be operated in two different modes: In the *mass spectrum mode* (MS mode) the chopper is moved completely out of the particle beam and the mass spectrum is scanned over a range of 1 to 300 amu with a frequency of 3 Hz. In this mode a quantitative measurement of the bulk mass concentrations of the volatile and semi volatile aerosol species is performed. For background measurement the particle beam is blocked by the chopper.

In the second mode, the *Time-of-Flight mode* (ToF mode), the particle beam is chopped by the chopper (2% duty cycle) and the time resolved MS signal is recorded for selected masses that are typical for some aerosol species. During the acceleration in the aerodynamic lens the particles reach a final velocity depending on their aerodynamic size: Larger particles are slower, smaller particles are faster. Using a particle velocity calibration (performed with laboratory particles of known sizes) the flight times, measured in the ToF mode can be converted into sizes, yielding to particle size distributions for the species that belong to the selected masses.

### **Field Operation**

During the PMTACS-NY 2001 summer intensive the AMS was operated at the measurement site at Queens College in Queens, New York from June 30<sup>th</sup> until August 5<sup>th</sup>. The measurement site was located at the edge of parking field # 6 adjacent to an athletic field.

The AMS was housed in an air-conditioned trailer together with other aerosol instruments. The sampling inlet was mounted on a tower at a height of 5 m above ground next to the trailer. The sample line from the inlet base to the instrument was approximately 2m. Ambient air was sampled via a PM<sub>2.5</sub> cyclone at 10 l/min and through 14.1 mm ID copper tubing. The tube diameter was chosen to minimize losses by impaction and gravitational settling for the given flow rate. The overall transport losses were calculated using simple formulas for diffusion, gravitational settling and impaction losses (Hinds (1999)) for the given inlet geometry. For 20 nm particles the losses are below 3 %, dominated by diffusion; at 2.5 µm the impaction and settling dominated losses are about 5 %. The minimum of losses is less than 0.3 % for 260 nm particles.

At the entrance of the AMS inlet, a sample flow of 0.4 l/min is isokinetically extracted from the 10 l/min transport flow. 0.3 l/min of this flow is diverted to a CPC (TSI 3025), while 0.1 l/min is introduced into the AMS inlet. The transport losses within the first few centimeters of the AMS inlet entrance are much larger than in the copper tube, especially for the small particles, due to the very low flow rate: 7 % for 20 nm particles and about 3 % for 2.5  $\mu$ m particles, with a minimum of 0.4 % around 450 nm. The range of the total inlet losses are about 10 % for 20 nm and 2.5  $\mu$ m particles, with a minimum of 0.7 % for 350 nm particles.

In the field the AMS was operated to periodically switch (i.e. 20s cycle) between the ToF mode (size information for pre-selected masses) and the MS mode (complete mass spectrum for bulk aerosol). In the ToF mode the size distribution information was measured for mass 15 amu (ammonium), 30 and 46 amu (nitrate), 32 amu (air, nitrogen), 48, 64 and 80 amu (sulfate) and 55, 57, 69, 71 and 123 amu (organics). In MS mode the mass spectrum was scanned from 1 to 300 amu at a 3 Hz rate. After 5 seconds of analyzing the particle vapor the particle beam was blocked for another 5 seconds by moving the chopper wheel completely into the beam to measure the background signal. Every 10 minutes the mass spectrum and size distribution averages were saved to disc.

The multiplier signal and gain, the inlet flow rate, and the efficiency of ionization and transmission of the ions through the QMS are necessary to calculate particle mass concentration. The inlet flow meter was calibrated in the lab before the campaign, while ionization and transmission efficiency calibrations were performed several times in the field. The calibration procedure involves producing ammonium nitrate particles of 350 nm mobility diameter and introducing them into the AMS. The average number of ions per particle reaching the multiplier entrance is determined from the multiplier signal, using the freshly calibrated multiplier gain. The ionization and transmission efficiency is the ratio of ions per particle as determined in this calibration procedure and molecules per particle, calculated from particle diameter, particle density and molecular weight. The multiplier gain was calibrated using single ions, produced from the background gas in the ionizer region at reduced electron current.

# **Data Processing**

Since the ionization and transmission efficiency didn't change significantly over the time of deployment, an average value was used for the whole campaign.

The mass concentration of a species is calculated from the multiplier signal, using the multiplier gain calibration and the ionization and transmission efficiency calibration together with the inlet flow meter reading.

For *nitrate* particles the mass concentration is simply the sum of the mass concentrations calculated for all fractions: 30 amu (NO) and 46 amu (NO<sub>2</sub>). Losses due to incomplete focusing of internally mixed nitrate and sulfate particles were accounted for by multiplication of the nitrate mass concentrations with a factor of 2.34. This factor was derived from the overall sulfate correction factor (2.9, determined from AMS – PILS comparison, see below) using the relative ionization efficiency of sulfate as determined in laboratory experiments.

Chloride typically is measured as Cl<sup>+</sup> and HCl<sup>+</sup>, considering chloride's two isotopes this includes mass 35, 36, 37 and 38 amu. As a result of strong interferences with fractions of organic species, only mass 35 and 36 amu are used to calculate the chloride mass concentration. The chloride mass concentration is corrected, using the natural isotope ratio of <sup>35</sup>Cl and <sup>37</sup>Cl.

Ammonium has fractions at mass 15, 16 and 17 amu (NH, NH<sub>2</sub>, NH<sub>3</sub>). Due to interferences with  $O_2^{2+}$  (16 amu) and fractions of water (OH, 17 amu) only the signal at 15 amu is used to calculate the ammonium signal. The omitted fractions are accounted for by multiplying the signal of mass 15 amu by a factor of 10, which was determined by the manufacturer by comparison with other instruments (PILS-IC).

Sulfate fractions (and their isotopes) are located at mass 48, 64, 65, 80, 81, 82 and 98 amu. Interferences with organic fractions at mass 65, 81, 82 and 98 amu requires that only mass 48, 64 and 80 amu be used to calculate the sulfate mass concentration. This mass concentration is multiplied by an empirical factor of 1.22 to account for the omitted masses. In addition sulfate particles exhibited a systematic loss of signal that thought to be related to incomplete focusing of the sulfate particles in the aerodynamic lens and to a less effective ionization of the sulfate particles in the ionizer. This loss in sulfate signal is accounted for by multiplying all sulfate mass concentrations with an empirical factor of 2.9. This factor was determined once by comparison of the AMS sulfate mass concentration with the sulfate mass concentration measured by another instrument (PILS).

Lab experiments by Aerodyne Research Inc. indicate that the ionization efficiency of organic molecules is larger than for inorganic molecules of equal molecular weight. This is accounted for by multiplication of all organics signals by a factor of 0.7. The *total organics* signal is calculated adding up all masses of the mass spectrum larger than 11 amu besides the masses that belong to fractions of other species or air.

The *size distribution data*, measured at the different masses were averaged for sulfate, nitrate, ammonium and total organics. The size bins of the raw data are determined by the time steps of the time resolved data acquisition. To obtain more useful size bins, that reflect the size resolution of the instrument, the raw data were re-binned, using 20 size bins per decade, equally spaced in the  $\log D_p$  space.

The size distribution raw data were transformed into mass concentration units ( $\mu g/m^3 \mu m$ ) and corrected to match the mass concentrations, measured by the MS mode.

# **Data Quality**

All valid data that were not qualified in any way were flagged with 'V0'.

Valid 10 min data that were below the detection limit were flagged with 'V1'. The detection limit was calculated for the MS mode using the noise of the background signal. For the size distributions the detection limit was calculated for every bin separately using the noise of the signal at large flight times (corresponding to large particle sizes above 2.5 µm, i.e. above the cut-off of the inlet cyclone) and the number of bins of the raw size distribution, accumulated in the re-binned size distribution bins. Accordingly every bin of every size distribution data set has its own flag.

Valid 10-min data, where the actual sampling and averaging time is below 7.5 min, but not below 5 min were flagged with 'V2'.

The data of one cycle were flagged with 'V6'. During this sampling cycle the heater was turned off.

Values missing, because the instrument was not sampling were flagged with 'M1'.

Values invalidated, because the sampling and averaging time was below 5 min were flagged with 'M2'.

No other flags were used.

### **Reference**

Jayne, J.T., Leard, D.C., Zhang, X., Davidovits, P., Smith, K.A., Kolb, C.E., Worsnop, D.R.: Development of an Aerosol Mass Spectrometer for Size and Composition Analysis of Submicron Particles. Aerosol Science and Technology **33** (2000) 49-70